

PARAMETERS AND TESTING TIME SPAN FOR HOMOEOPATHIC DRUGS

Drug Name	Parameters	Time Span
MOTHER TINCTURE	<ul style="list-style-type: none"> • pH • wt/ml • Total Solid • Alcohol Content • Identification of Plant/Animal/Chemical • Absorption maxima 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days</div>
DILUTION	<ul style="list-style-type: none"> • Alcohol Content • Test for Aldehyde, Ketone, Fusel oil, Methanol 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days</div>
TABLET	<ul style="list-style-type: none"> • Ash value • Lactose test • Disintegration • Identification of drug • Other test, if any 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days</div>
FORMULATION	<ul style="list-style-type: none"> • pH • Alcohol Content • Identification of drugs • Sterility Test 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days + 7-10 days for Microbial analysis</div>
EYE DROP	<ul style="list-style-type: none"> • pH • Alcohol Content • Identification • Sterility Test 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days + 7-10 days for Microbial analysis</div>
EAR DROP	<ul style="list-style-type: none"> • pH • Alcohol Content • Identification 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days</div>
GLOBULES	<ul style="list-style-type: none"> • Size • Porosity • Sugar content • Solubility • Ash value • Foreign matter 	<div style="border: 1px solid black; background-color: #cccccc; padding: 5px; width: fit-content; margin: 0 auto;">2 working days</div>

OINTMENT	<ul style="list-style-type: none">• Homgenity• Clarity• Identification of drug	2 working days
OIL	<ul style="list-style-type: none">• Identification of drug• Determination of Acid value, Saponification value, Iodine value	2 working days